

# iWash™ Slide Cleaner

“The innovative iWash Slide Cleaner has integrated into our cell culture workflow seamlessly, it is easy to use and adopted well by researchers. There is minimal maintenance, uses minimal bench space and helps towards our sustainable targets. The cost benefits by reusing slides allows our funding to go further, as well as being eco-friendly.”

**Vipul Bhakta**  
Queens University  
London, UK

## For automated cell counting slides

*TekMatic offers an easy-to-use slide cleaner as a way of minimizing the cost of automated cell counting and reducing biohazardous waste. Using patent pending direct injection technology, the iWash Slide Cleaner completely cleans and dries a slide in less than 30 seconds. The iWash Slide Cleaner is compatible with a variety of popular image-based cell counters.*



## Compatibility

- Bio-Rad TC20™
- Invitrogen Countess II™
- NanoEnTek Eve™ & Arthur™
- Logos Biosystems Luna™ Series
- Olympus R1™
- Nexelom Cellometer™ Series

## Features

- Wash cell counter slides in less than 30 seconds
- Reuse the same disposable cell counter slide up to 20 times
- Completely automated
- Reduce biohazardous waste
- ROI in less than 3 months

