


# ODTC Versions

## XL Version and Compact Version

The 96 XL gives you plate flexibility, you can use the high volume PCR plates and also the standard volume PCR plates.  
 The 96 XL & 384 XL provide more gripper space, but the gripper space of the compact ODTC is usually sufficient.  
 The compact ODTC leaves even more space for other devices on your liquid handling workstation.


### XL gripper space





### Ventilation Outlet Versions:

ventilation outlet to the back, to the left, to the right, and downwards


**all options on deck**  
**freely chose the ideal position**




**air exhaust downwards**



ODTC XL, 113mm



ODTC Compact, 103mm



## ODTC Scope of Supply

all items required to operate the ODTC are included in the scope of supply.

included are Power & Control Unit, PMS Software, Script Editor, Sealing Cover, Cables & other accessories.

## ODTC Prices

Price Category A are all identical, Price Category B are all identical

## ODTC Versions

Thermal Cycler						
Product	INHECO #	Plates	Lenght	Width	Height	Price Category
<b>ODTC® 96 XL Left</b> (ventilation left) for all 96 well PCR plates (150µL and 300µL)	8100300	300µL 150µL	285.2	156.5	129.3	A
<b>ODTC® 96 XL Back</b> (ventilation back) for all 96 well PCR plates (150µL and 300µL)	8100301	300µL 150µL	285.2	156.5	129.3	A
<b>ODTC® 96 XL Right</b> (ventilation right) for all 96 well PCR plates (150µL and 300µL)	8100302	300µL 150µL	285.2	156.5	129.3	A
<b>ODTC® 96 XL Down</b> (ventilation down) for all 96 well PCR plates (150µL and 300µL)	8100303	300µL 150µL	285.2	156.5	129.3	A
<b>ODTC® 96 Left</b> (ventilation left) only for max 150µL well PCR plates	8100100	150µL	248	156.5	124.3	A
<b>ODTC® 96 Back</b> (ventilation back) only for max 150µL well PCR plates	8100101	150µL	248	156.5	124.3	A
<b>ODTC® 96 Right</b> (ventilation right) only for max 150µL well PCR plates	8100102	150µL	248	156.5	124.3	A
<b>ODTC® 96 Down</b> (ventilation down) only for max 150µL well PCR plates	8100103	150µL	248	156.5	124.3	A
<b>ODTC® 384 XL Left</b> (ventilation left) XL for extended gripper space	8100400	384 well	285.2	156.5	129.3	B
<b>ODTC® 384 XL Back</b> (ventilation back) XL for extended gripper space	8100401	384 well	285.2	156.5	129.3	B
<b>ODTC® 384 XL Right</b> (ventilation right) XL for extended gripper space	8100402	384 well	285.2	156.5	129.3	B
<b>ODTC® 384 XL Down</b> (ventilation down) XL for extended gripper space	8100403	384 well	285.2	156.5	129.3	B
<b>ODTC® 384 Left</b> (ventilation left)	8100200	384 well	248	156.5	124.3	B
<b>ODTC® 384 Back</b> (ventilation back)	8100201	384 well	248	156.5	124.3	B
<b>ODTC® 384 Right</b> (ventilation right)	8100202	384 well	248	156.5	124.3	B
<b>ODTC® 384 Down</b> (ventilation down)	8100203	384 well	248	156.5	124.3	B